

GENERAL LABORATORY RISK ASSESSMENT

1. General information.

| | | | |
|--|-----------------------------------|------------------|---------------------|
| Department: | Department of Nanotechnology | Date: | October 23, 2024 |
| Procedure: | Lp-PLA2 ECL Biosensor Fabrication | Revision Number: | 1 |
| Principal Investigator: | Dr. John Saad | PI Phone Number: | 519-888-4567 x39296 |
| PI Signature: *By signing you are indicating that the tasks are planned for in such a way that the risk is tolerable. <div style="text-align: right; margin-top: 20px;"> <u>X</u> </div> | | | |

2. Identify if any of the following hazards or materials are present.

| | | | |
|---|--|---|---|
| <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | Nanomaterials | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | Cannabis |
| <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | Biohazards | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | Nonionizing Radiation |
| <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | Radioactive materials | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | Designated Substances (acrylonitrile, benzene, silica, isocyanates, vinyl chloride, As, Pb, Hg, etc.) |
| <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | X-ray sources | | |
| <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | Class 3B or Class 4 Lasers | | |

If you checked “yes” for any of the items above, review the associated program to ensure you have managed those requirements prior to or in conjunction with completing this risk assessment.

3. Describe the project steps in point form detail. Identify equipment and materials at relevant steps.

A) CoFe PBA Fabrication

1. Solution A is made by dissolving 436.2 mg of cobalt (II) nitrate hexahydrate ($\text{Co}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$) and 441.2 mg of trisodium citrate dihydrate ($\text{C}_6\text{H}_5\text{Na}_3\text{O}_7 \cdot 2\text{H}_2\text{O}$) in 50 mL of deionized water
2. Solution B is made by dissolving 329.25 mg of potassium ferricyanide (III) ($\text{K}_3[\text{Fe}(\text{CN})_6]$) in 50 mL of deionized water
3. The two solutions were mixed with a magnetic stir plate for 5 minutes
4. The new solution was aged for 24 hours at room temperature
5. The solution was spun in a centrifuge at 8000 rpm for 5 minutes
6. The dark purple precipitates were collected and washed three times with deionized water and ethanol separately
7. The substance was dried overnight (12 hours) at 60°C

B) PEI@CoFe PBA Combination

1. The CoFe PBA (10 mg) were dispersed via sonification in 10 mL of deionized water alongside 100 μL of polyethylenimine (PEI)
2. The mixture was magnetically stirred for 4 hours at 800 rpm at room temperature
3. The mixture was spun in a centrifuge at 8000 rpm for 5 minutes to separate the precipitates
4. The precipitates were washed three times with deionized water
5. The precipitates were dispersed in about 10 mL of deionized water (subject to yield) to obtain 1 mg/mL of PEI@CoFe PBA
6. The mixture was stored at 4°C in the dark

C) AuNP@CoFe PBA Nanocomposites Adhesion

1. 1 mL of the PEI@CoFe PBA solution was added to 9 mL of a 25 nm AuNP colloid and stirred at room temperature for 2 hours
2. The solution was spun in a centrifuge at 8000 rpm for 5 minutes
3. The precipitate was collected and washed three times with deionized water
4. The particles were dispersed in about 5 mL of deionized water (subject to yield) to make 0.2 mg/mL of AuNP@CoFe PBA nanocomposites

D) ITO Electrode Preparation

1. The ITO-coated substrate was sanitized before usage

- a. For ITO-coated glass/PEN, it was sonicated in a 1:1 v/v 2% sodium hydroxide (NaOH) and ethanol for 30 minutes; acetone and deionized water were also used for sonification for 30 minutes
2. It was sonicated in a 30% ammonia solution for 30 minutes and soaked overnight to hydrolyze its surface
3. The ITO was cleaned with deionized water and dried with nitrogen gas 3 times
4. 10 μL of a 0.01% v/v APTMS and ethanol solution was dropped onto the ITO surface
5. It was placed in a water saturated environment at 55°C for 4 hours
6. 50 μL of the AuNP@CoFe PBA solution was dropped onto the ITO surface for 5 minutes
7. The excess solution was washed off with deionized water and dried with nitrogen gas
8. The electrode was stored at 4°C until next stage

E) Antibody Adherence

1. The 10 μL Lp-PLA2 antibody mixture (0.2 ng/mL) was dropped onto the electrode surface and incubated at 30°C for 2 hours
2. The excess solution was washed off with 0.01 M of phosphate buffered saline (PBS) at 7.4 pH
3. The remaining exposed active sites were closed off with 10 μL of 1% v/v bovine saline albumin (BSA) at 25°C for 1.5 hours
4. The excess solution was washed off with 0.01 M of PBS at 7.4 pH and dried with nitrogen gas
5. The completed sensor was stored at 4°C until usage

F) Sample Testing

1. A 0.2 M PBS solution at 8 pH containing 1.0×10^{-5} M of Luminol was prepared as the electrolytic solution
2. The ITO electrode was placed on the circuit containing the counter and reference electrodes
3. The electrolytic solution was dropped onto the circuit ensuring all the electrodes are covered
4. The circuit was placed into the dark box containing the photodiode sensor
5. The circuit was powered by a pulse square wave with a V_{pp} of 1.7 V with a maximum of 1.3 V and minimum of -0.4 V and a period of 3 seconds
6. Its baseline light intensity was recorded, and solution was dried off with nitrogen gas
7. The sample containing the Lp-PLA2 was coated on the ITO sensor and left to incubate at 35°C for 80 minutes maximum
8. The electrolytic solution was dropped onto the circuit ensuring all the electrodes are covered
9. The circuit was placed into the dark box containing the photodiode sensor
10. The circuit was powered by a pulse square wave with a V_{pp} of 1.7 V with a maximum of 1.3 V and minimum of -0.4 V and a period of 3 seconds
11. Its light intensity was recorded

4. Identify the WHMIS hazard classes, categories, and anticipated control measures used to reduce worker exposure. (Hover mouse here to learn how to add more rows)

| Name of chemical | List all WHMIS physical and health hazard classes and categories (drop-down list and free-form text) | Category / Type | Using WHMIS precautionary statements as a guide, identify what control practices are required to minimize worker exposure for handling these chemicals – consider ventilation, PPE, containment, etc. |
|---------------------------------------|--|-----------------|---|
| 3-Aminopropyltrimethoxysilane (APTMS) | Skin corrosion | 1 | <ul style="list-style-type: none"> Chemical must ALWAYS be used in the fume hood or a well-ventilated environment Chemical can react with water and produce <u>methanol vapours</u> (BAD) so use in fume hood ALWAYS wear eye protection and gloves around substance |
| | Skin irritation | 2 | |
| | Serious eye damage | 1 | |
| | Eye irritation | 2 | |
| | Substances and mixtures which, in contact with water, emit flammable gases | 1 | |
| Cobalt (II) nitrate hexahydrate | Oxidizing solids | 2 | <ul style="list-style-type: none"> Keep out of reach of flammable material and chemicals when possible Chemical must ALWAYS be used in the fume hood or a well-ventilated environment ALWAYS wear eye protection and gloves around substance <ul style="list-style-type: none"> ENSURE gloves are the following materials: <ul style="list-style-type: none"> Natural Rubber Nitrile Rubber Neoprene PVC |
| | Acute toxicity - oral | 4 | |
| | Skin sensitizer | 1 | |
| | Serious eye damage | 1 | |
| | Acute toxicity - inhalation | 4 | |
| | Respiratory sensitizer | 1 | |
| | Germ cell mutagenicity | 2 | |
| | Carcinogenicity | 1B | |
| | Reproductive toxicity | 1B | |
| | Skin corrosion | 1B | |
| | Specific target organ toxicity - repeated exposure (list organs) <ul style="list-style-type: none"> Cardiovascular Organs Thyroid Bone Marrow | 1 | |
| Potassium ferricyanide (III) | Acute toxicity - oral | 4 | <ul style="list-style-type: none"> Will decompose into <u>hydrogen cyanide</u> (BAD) when heated, so NEVER heat this outside a fume hood or well-ventilated area Will decompose into <u>cyanide</u> (ALSO BAD) when reacted with an acid, so NEVER bring into contact with an acid Chemical must ALWAYS be used in the fume hood or a well-ventilated environment ALWAYS wear eye protection and gloves around substance |
| | Skin irritation | 2 | |
| | Eye irritation | 2 | |
| | Specific target organ toxicity - single exposure (list organs) <ul style="list-style-type: none"> Respiratory Organs | 3 | |
| | Reproductive toxicity | 2 | |
| Polyethylenimine (PEI) | Acute toxicity - oral | 4 | <ul style="list-style-type: none"> ALWAYS wear eye protection and gloves around substance |
| | Serious eye damage | 1 | |
| Indium tin oxide | Specific target organ toxicity - repeated exposure (list organs) <ul style="list-style-type: none"> Respiratory Organs | 1 | <ul style="list-style-type: none"> Do not inhale ITO dust, so minimize chances of producing this Wash hands thoroughly after handling material |
| | | | |
| Bovine saline albumin (BSA) | Specific target organ toxicity - single exposure (list organs) <ul style="list-style-type: none"> Respiratory Organs | 3 | <ul style="list-style-type: none"> Do not inhale fumes or vapours from substance, so handle in well-ventilated area |
| Luminol | Acute toxicity - oral | 4 | <ul style="list-style-type: none"> Wear eye protection and gloves around substance Handle in well-ventilated area |
| | Skin irritation | 2 | |
| | Eye irritation | 2 | |
| | Specific target organ toxicity - single exposure (list organs) <ul style="list-style-type: none"> Respiratory Organs | 3 | |

| | | | |
|------------------|---|----|--|
| Sodium hydroxide | Corrosive to metals | 1 | <ul style="list-style-type: none"> Produces heat when dissolved in water or neutralized with acid, so be wary when using with flammable substances ALWAYS wear eye protection and gloves around substance |
| | Skin corrosion | 1A | |
| | Serious eye damage | 1 | |
| | Eye irritation | 2 | |
| Ethanol | Flammable liquids | 2 | <ul style="list-style-type: none"> Be cautious of placement near an ignition source or an oxidizer Wear eye protection around this substance |
| | Eye irritation | 2 | |
| Acetone | Flammable liquids | 2 | <ul style="list-style-type: none"> Be cautious of placement near an ignition source or an oxidizer Work with substance in a well-ventilated environment Wear eye protection around this substance |
| | Eye irritation | 2 | |
| | Specific target organ toxicity - single exposure (list organs) <ul style="list-style-type: none"> Narcotic Effects | 3 | |
| Ammonia solution | Skin corrosion | 1B | <ul style="list-style-type: none"> ALWAYS wear eye protection and gloves around substance Work with substance in a well-ventilated environment |
| | Specific target organ toxicity - single exposure (list organs) <ul style="list-style-type: none"> Respiratory Organs | 3 | |

5. List equipment being used.

Equipment used consists of devices that impart energy or contain reactions – examples include rotovaps, ovens, pressure devices, material test stands etc. ([Hover mouse here to learn how to add more rows](#))

| Identify equipment (Click box to enter text) | Step # | List Potential Equipment Hazards | Identify anticipated controls to control the identified risks |
|---|------------------|---|--|
| Magnetic Stir & Hot Plate | A3, A7, B2, C1 | Vibration | <ul style="list-style-type: none"> Place a sign nearby the hot plate to signify its usage |
| | | Source of heat - potential burns | |
| | | Moving parts | |
| Sonicator | B1, D1, D2 | Vibration | <ul style="list-style-type: none"> Use ear protection due to high frequency noise generated Place a sign nearby the sonicator to signify its usage |
| | | Noise | |
| Centrifuge | A5, B3, C2 | Vibration | <ul style="list-style-type: none"> ALWAYS ensure the inner stage is balanced with 2 or more vials Place a sign nearby the centrifuge to signify its usage |
| | | Moving parts | |
| Function Generator | F5, F10 | Unprotected electrical source | <ul style="list-style-type: none"> Ensure it is not set to a high voltage before touching outputs |
| UV-Vis | Characterization | Radiation emission - UV, IR, Visible light, microwave, radiofrequency | <ul style="list-style-type: none"> Do not look at light source during operation |
| SEM | Characterization | Moving parts | <ul style="list-style-type: none"> Follow approved operation steps to be safe |
| | | High pressure release | |

6. Hazardous conditions.

Please identify how experimental conditions (pressure, temperature, humidity, etc.) may alter the behaviour of chemicals being used. Remember to consider reactive intermediates as well.

- APTMS will be dropped onto the ITO electrode and placed in a water saturated environment at 55°C for 4 hours, which will produce methanol as a byproduct
 - This is very dangerous, so it must be done in the fume hood and the experimenter must wear proper PPE
- Residual potassium ferricyanide (III) may be present during the overnight drying process at 60°C, which may decompose and produce hydrogen cyanide gas
 - This is very dangerous, so this step must also be left in the fume hood with a corresponding sign and the experimenter must be careful not to inhale the gas

7. Managing hazardous waste.

Please include waste disposal methods in your SOPs. Guidance can be found on the [Hazardous Waste website](#) or by emailing esf@uwaterloo.ca. (Hover mouse here to learn how to add more rows)

| Contents and Anticipated Class | Waste Type | Anticipated Amount |
|---|--|--------------------|
| Potassium nitrate (Aqueous liquid waste – inorganic salts) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | <100 mL |
| Polyethylenimine (Aqueous liquid waste – organics) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | <10 mL |
| Sodium hydroxide (Base) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | 250 mL |
| Ethanol (Organic solvents) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | 400 mL |
| Acetone (Organic solvents) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | 500 mL |
| Ammonia (Aqueous liquid waste – inorganic salts) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | 500 mL |
| Phosphate buffered saline (Aqueous liquid waste – inorganic salts) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | 300 mL |
| Luminol (Organic solvent) | <input type="checkbox"/> Solid <input type="checkbox"/> Radioactive <input type="checkbox"/> Biological <input checked="" type="checkbox"/> Liquid <input type="checkbox"/> Battery <input type="checkbox"/> Chemical | 50 mL |
| <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | All necessary labels, containers, transportation means are available to start the research process. | |
| <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | All waste generators have taken the online Chemical Waste Segregation SO2070 | |

8. Standard operating procedures and emergency planning.

All medium to high-risk activities require an SOP. Work with toxic, pyrophoric, or water reactive materials require emergency planning SOPs. Identify what SOPs will be created for this project in the table below and where they are located. [SOP template is available from the Safety Office](#).

| SOP Name | Procedure available | Indicate how this SOP covers anticipated risks |
|--------------------------------------|---|--|
| Overall Process | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <ul style="list-style-type: none"> Guides the experimenter safely Ensures nothing is done incorrectly |
| Spill, Exposure, or Emergency | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <ul style="list-style-type: none"> In the scenario that something occurs, the experimenter can safely navigate the situation Mainly deals with spillage, and accidental exposure |

9. Personal protective equipment.

Note: Closed toed shoes and lab coat are mandatory for work with chemicals.

| PPE Type | PPE Storage Location | When it is Worn |
|----------------|------------------------|-----------------|
| Nitrile Gloves | QNC Undergraduate Labs | Always |
| Goggles | Personal Item | Always |

10. Supervisor commentary.

11. Worker sign-off.

By signing the sheet below, you acknowledge that you have:

1. Understood the stipulations, hazardous, and control requirements outlined in this document.
2. You have completed practical training and had the opportunity to ask questions

| Name (Print) | Signature | Date |
|----------------------|-----------|-------------------------------|
| Brandon Kong | X | Click or tap to enter a date. |
| Magilan Varatharuban | X | Click or tap to enter a date. |
| Rand Dakhil | X | Click or tap to enter a date. |
| Andrew Wang | X | Click or tap to enter a date. |